

Product Description

Procollagen Dressing is an extracellular matrix (ECM) dressing that is available from Day One. It retains the structure and function of the intact, native ECM4 to help support the patient's degraded matrix. It also demonstrates broad-spectrum MMP activity reduction, especially against critical collagenases, gelatinases, and additional MMPs.

Features

- **Available from Day One** - Procollagen Dressing is available from Day One. It retains the structure and function of the intact, native ECM4 to help support the patient's degraded matrix. It also demonstrates broad-spectrum MMP activity reduction, especially against critical collagenases, gelatinases, and additional MMPs.

ENDOTERMAL DERMAL TEMPLATE
Procollagen Dressing is available from Day One. It retains the structure and function of the intact, native ECM4 to help support the patient's degraded matrix. It also demonstrates broad-spectrum MMP activity reduction, especially against critical collagenases, gelatinases, and additional MMPs.

It also demonstrates broad-spectrum MMP activity reduction, especially against critical collagenases, gelatinases, and additional MMPs.